

# RVA world

The technical journal of Newport Scientific

JANUARY 2001

ISSUE 17

## Use of the RVA™ in the Evaluation of Pasta Processing

Paul J. Whalen

Whalen Consulting Inc.,  
200 Fifth Street NW, Ste. L,  
Elk River, MN 55330, USA




Pasta is a half-product intended for further preparation at home or by food service. In this sense it is similar to other products such as ready-to-eat breakfast cereals and snacks wherein the half-product or pellet is then puffed, flaked or fried. Thus, the ingredient and process effects should be reflected in the half-product qualities. This work was initiated by observations made from RVA analysis of pasta product purchased off the shelf at retail grocery stores. Standard RVA profiles (25°C-95°C-25°C) showed a product profile (Figure 1) reflecting low-cook characteristics with essentially no difference in peak paste between domestic and imported pasta brands. However, application of a critical paste viscosity (Figure 2) showed extremely different profiles and speculation as to the cause — the durum wheat, semolina, extrusion or drying regimen.

This project was initiated to look for the cause of the difference in C-paste RVA. The RVA was used to observe the differences in the ingredients and the effects of extrusion, as well as the effects of the different possible drying regimens. The RVA methods have been previously published (P.J. Whalen, 1998. *Cereal Foods World*, Vol. 43, No.2:69-72).

### Figure 1

Standard RVA profiles of domestic and imported pasta brands.


### RVA™ PASTA PRODUCTS — STD PROFILE



### Figure 2

Critical paste RVA profile of imported and domestic pasta (spaghetti).

### RVA™ DIFFERENTIATION OF PASTA PRODUCTS



## RVA™ Methods

All samples were run on a dry weight basis (dsb), adjusting for moisture.

### Standard RVA Profile

Foss Std1 method extended to 15 min.

Idle = 50°C; mix 10 sec. @ 960 rpm; run at 160 rpm;  
Temp. = hold 50°C, 1.0 min.;  
Ramp to 95°C at 4.42 min.;  
Hold 95°C at 7.12 min.;  
Cool at 11 min. to 50°C;  
End at 15 min.

### RVA Sweep

Idle = 50°C; mix at 960 rpm, 10 sec.; run at 160 rpm;  
Ramp to 95°C in 15 min. (End)

### Critical Paste RVA

Idle = 25°C;  
Mix at 960 rpm, 10 sec.; run at 160 rpm;  
Ramp to 63°C at 5 min.;  
Hold 63°C, Ending at 15 min.

## Wheat

The wheat was purchased for this study through North Dakota State University (NDSU), Fargo, North Dakota, USA. The strong and weak gluten durum wheat was evaluated by RVA. No durum wheat sample was available for the commercial semolina. All moistures were comparable (13.8 - 14%). Mixograph values for the semolina were: strong = 7; commercial = 4; and weak = 1.

### IN THIS ISSUE...

- Recent references
- Method lift-out
- Tips & tricks

**Std profile — RVA Data**  
(in centipoise)


Wheat	Peak	Trough	Setback*
Strong	1002	599	1427
Weak	1535	1026	2323

\*Setback as referred to herein is the final viscosity.

As shown in Figures 3a and 3b below, the strong gluten wheat was differentiated by both the standard RVA profile and the sweep method. The standard profile showed very large differences at all key points of the profile. The weak gluten durum was 50% higher in peak paste, 71% higher at the trough, and had a 63% higher setback.


**Figure 3a** Durum wheat RVA profiles of strong and weak gluten durum via sweep profile.

**DURUM WHEAT — SWEEP METHOD, STRONG & WEAK**



**Figure 3b** RVA sweep profiles for the semolina products.

**RVA™ SWEEP PROFILES — SEMOLINA: STRONG, WEAK & COMMERCIAL**



**Dryer Profiles**

The profiles used for the low temp. (40°C, 19 - 20 h), high temp. (72°C, 12 h) and ultra-high temp. (89°C, ~7 h) drying were the standard conditions (humidity profile) used at the NDSU Durum/Pasta Lab. Samples were taken exiting the extruder, then just prior to starting the dryer (zero), and at key transitions in accordance with each profile (low = 1 h, 11 h, and final; high = 4 h, 10 h, and final; ultra-high temp. = 1 h and final). Standard RVA profiles were run on all samples and sweep profiles made on the final products to identify a critical paste temperature.


**Results**

**Effects of Extrusion**

An extruded product will generally always reflect a degradative RVA profile compared to the starting material. RVA profiles can follow this effect very well since the profile decreases and shifts as the material is changed by mechanical input (Figure 4). This is essentially what happened in the pasta process along with some other significant effects.

**Figure 4** Effects of extrusion on strong gluten semolina — initial product transformation.

**EXTRUSION EFFECT — STRONG GLUTEN, SEMOLINA vs EXTRUDER**



**Effect of Drying**


The trend during the drying regimen is for the peak paste to shift left as the temperature increases during drying as shown for the high temp. (HT) treatment in Figure 5. This is evident in the ultra-high temp. (UHT) drying profile as well. Although only one sample could be taken after starting the dryer, it is evident that the

initial shift is quite dramatic. Also, the UHT setback was always higher in the final pasta product for all semolina varieties. This further indicates cross-linking and similar reactions as seen in thermally pre-conditioned materials. These may or may not be associated with protein interactions. It is assumed and accepted as general knowledge that the high relative humidity aids in the conductive heat transfer and assists in the drying effect as in other half-product treatments. Figure 5 shows the UHT profiles for the strong gluten. The commercial source behaved very similar to the strong gluten.

The C-paste method confirms and verifies that the differences observed in the off-the-shelf samples of pasta are due to the different drying profiles employed — low temp. vs ultra-high temp. In Figures 6 and 7, the pattern shown for the different sources of durum semolina by the C-paste method indicate that the material is continuing to cook in the dryer at the UHT temperatures. This is why the absorption changes are so radical and are so easily detected by RVA C-paste for UHT pasta. UHT drying definitely continues the gelation of starch and this in turn would promote an interaction/network with any protein complex. The fact that the material retrogrades so readily is also indicative of the structure being formed in extrusion and drying. That Pollini (Pasta and Noodle Technol., AACC, 1996) notes a broader performance capability and tolerance for poorer grades of semolina by this method of drying is indicative of the interplay between starch cooking in the dryer and the extrusion forming process.

**Figure 5** Trend in peak paste shift at high temp. drying regimen for strong gluten pasta.


**TREND IN PEAK PASTE — STD RVA™, STRONG GLUTEN, HT**



A comparison of the general effects of the drying regimen is shown in Figures 6 and 7. UHT showed the highest paste viscosity followed by the high temp., and the lowest being the low temp. profile. These differences are directly related to the drying conditions (especially temperature) and the structure set-up in the pasta by the dryer. This trend was found to be the same for all three semolinas with varying degrees of differences among semolinas.


**Figure 6**  
Critical paste profile for the strong gluten semolina (63°C profile).

**CRITICAL PASTE RVA™**  
— STRONG GLUTEN PASTA



**Figure 7**  
Critical paste profile for the weak gluten and commercial semolina (63°C).

**CRITICAL PASTE RVA™**  
— WEAK GLUTEN PASTA



## Conclusions

### Durum Wheat

Differences between the durum (strong and weak gluten) are very apparent via Std RVA.

### Semolina

Differences are apparent via both the standard RVA and the RVA sweep profile.

### Pasta

*Extruder* — standard RVA showed differences between semolinas;

*Dryer* — differences were found between semolinas and dryer profiles;

*Final Product* — differences were strongly differentiated via dryer temperature profile (low, high, UHT) using the critical paste RVA and by semolina type.

This work also points to mechanisms occurring in the pasta process that are important and include the different responses of the strong and weak gluten and the totally different response of the commercial grade semolina. There is no doubt about protein quality as a factor but the interaction and effect of starch is also important. The very reason the UHT does what Pollini indicates with lower quality semolina is starch based. It is entirely likely that there is an interactive mechanism at work but the UHT is driven by starch gelation in the dryer as a cook effect. Translation to sensory and other product quality factors needs to be shown. However, both the extruder and dryer are driving factors for manipulating sensory responses along with semolina quality. In this work, we have only taken the 'known' and monitored it. The real opportunity is to take what was uncovered and apply this information in a knowledgeable manner.

## TIPS & TRICKS

### Exporting Experiment Data

A question which is frequently asked is how to send data, results and graphs from our Thermocline for Windows (TCW) software to other Windows programs such as MS Excel or Word. Using TCW version 2.0 or later, this is easily achieved:

First, run TCW and press F9 to cross to the Collection Tool. Now go to the Options menu and check the option to 'Save Data Comma Delimited'.

TCW will now save two data files for every test that you run. The first is the normal dat file (e.g. mydata.dat) for display in TCW's Analysis Tool. The second is a text copy of your data with the extension csv (e.g. mydata.csv), in Comma Separated Variable format.

The csv file contains date, instrument, calibration and profile information, and the test data arranged in four columns under the heading 'Datapoints'. The data columns are time in seconds, viscosity in cP, temperature in 0.05°C increments and speed in RPM.

The csv file can be imported into other Windows programs.

For example, to load it into Excel 97, select File – Open – Files of type: Text Files – and enter the file name.

If you have already created TCW data files without corresponding csv files, you can create text files from these later by using the datatext.exe utility located in your TCW directory.

## REFERENCES

### Recent References Mentioning the RVA™

- AACC, 1998. General Pasting Method for Wheat or Rye Flour using the Rapid Visco Analyser. AACC Method 76-21, First Approval 10-15-97, Approved Methods of Analysis, 9th ed., Amer. Assoc. Cereal. Chem., St. Paul MN.
- Allan, A.M., Blakeney, A.B., Batten, G.D., and Dunn, T.S. 1999. Impact of grinder configurations on grinding rate, particle size, and trace element contamination of plant samples. *Commun. Soil Sci. Plant Anal.* 30(15&16):2123-2135.
- Allen, A., Blakeney, A., Batten, G., and Dunn, T. 1999. Grinding plant samples: Rate, recovery and particle size. In Proc. 49th Aust. Cereal Chem. Conf., eds J.F. Panozzo, M. Ratcliffe, M. Wootton and C.W. Wrigley, Royal Aust. Chem. Inst., Melbourne. (In press.)
- Allen, H.M., Craze, T.L., and Fleming, D.K. 1997. Selection of wheat suitable for steam buns. In Cereals '97, Proc. 47th Australian Cereal Chem. Conf., eds A.W. Tarr, A.S. Ross and C.W. Wrigley, pp. 87-90. Royal Aust. Chem. Inst., Melbourne.
- Batey, I.L., and Curtin, B.M. 1997. Through thick and thin — Factors affecting the viscosity of wheat starch. In Cereals '97, Proc. 47th Australian Cereal Chem. Conf., eds A.W. Tarr, A.S. Ross and C.W. Wrigley, pp. 369-372. Royal Aust. Chem. Inst., Melbourne.
- Corke, H. 1997. Manipulating starch quality for Asian food applications: Genetic, physical and chemical approaches. In Cereals '97, Proc. 47th Australian Cereal Chem. Conf., eds A.W. Tarr, A.S. Ross and C.W. Wrigley, pp. 103-105. Royal Aust. Chem. Inst., Melbourne.
- Crosbie, G.B., Ross, A.S., and Chiu, P.C. 1999. Starch and protein quality requirements of Japanese alkaline noodles (Ramen). *Cereal Chem.* 76(3):328-334.
- Crosbie, G.B., Ross, A.S., and Chiu, P.C. 1999. Effects of  $\alpha$ -amylase on flour paste viscosity measurements and relationships with alkaline noodle texture. *RVA World* 15:2-3.
- Dines, J., and Taylor, T. 1999. Frost, late maturity  $\alpha$ -amylase and Australian wheat receival standards. In Proc. 49th Aust. Cereal Chem. Conf., eds J.F. Panozzo, M. Ratcliffe, M. Wootton and C.W. Wrigley, Royal Aust. Chem. Inst., Melbourne. (In press.)
- Elliott, B. 1999. Análisis rápido de viscosidad. *Techno Food.* (1999):36-40.
- Fairbrother, A.H., and Birkett, J. 1997. Wheat quality associated with the extrusion of four wheat varieties and wheat flour/lupin flour mixes. In Cereals '97, Proc. 47th Australian Cereal Chem. Conf., eds A.W. Tarr, A.S. Ross and C.W. Wrigley, pp. 151-156. Royal Aust. Chem. Inst., Melbourne.
- Graybosch, R.A., Guo, G., and Shelton, D.R. 2000. Aberrant Falling Numbers of waxy wheats independent of  $\alpha$ -amylase activity. *Cereal Chem.* 77(1):1-3.
- Harris, D.J., Sipsas, S., and Petterson, D.S. 1997. Effect of milling procedures on recovery, protein and moisture content in grain legume flours. In Cereals '97, Proc. 47th Australian Cereal Chem. Conf., eds A.W. Tarr, A.S. Ross and C.W. Wrigley, pp. 350-354. Royal Aust. Chem. Inst., Melbourne.
- Jun, W.J., Seib, P.A., and Chung, O.K. 1998. Characteristics of noodle flours from Japan. *Cereal Chem.* 75(6):820-825.
- Kiribuchi-Otobe, C., Yanagisawa, T., Yamaguchi, I., and Yoshida, H. 1998. Wheat mutant with waxy starch showing stable hot paste viscosity. *Cereal Chem.* 75(5):671-672.
- Kruger, J.E. 1997. Wheat enzymes as a quality determinant. In Proc. Int. Wheat Quality Conf., eds J.L. Steele and O.K. Chung, pp. 101-106. Grain Industry Alliance, KS.
- Maningat, C.C., and Seib, P.A. 1997. Update on wheat starch and its uses. In Proc. Int. Wheat Quality Conf., eds J.L. Steele and O.K. Chung. Grain Industry Alliance, KS.
- O'Shea, M.G., Batey, I.L., Samuel, M.S., Curtin, B.M., Konik, C., and Morell, M.K. 1997. Analytical tools for determining natural variation in wheat starch structure and function. In Cereals '97, Proc. 47th Australian Cereal Chem. Conf., eds A.W. Tarr, A.S. Ross and C.W. Wrigley, pp. 143-147. Royal Aust. Chem. Inst., Melbourne.
- Pehm, S., Jackson, K.L., and Azudin, N. 1999. Evaluation of wheat quality parameters that optimize binding characteristics of shrimp feed pellets. In Proc. 49th Aust. Cereal Chem. Conf., eds J.F. Panozzo, M. Ratcliffe, M. Wootton and C.W. Wrigley, pp. 191-196. Royal Aust. Chem. Inst., Melbourne.
- Prakash, M., Ravi, R., and Rao, P.H. 1999. Effect of blending of raw and steamed wheat flour on the rheological characteristics of dough. *Eur. Food Res. Technol.* 210:119-122.
- Sasaki, T., Yasui, T. and Matsuki, J. 2000. Effect of amylose content on gelatinization, retrogradation, and pasting properties of starches from waxy and nonwaxy wheat and their F1 seeds. *Cereal Chem.* 77(1):58-63.
- Seib, P.A. 1997. Wheat starch as a quality determinant. In Proc. Int. Wheat Quality Conf., eds J.L. Steele and O.K. Chung, pp. 61-82. Grain Industry Alliance, KS.
- Sharma, R., Cornish, G.B., Jenner, C.F., Palmer, G.A., and Rathjen, A.J. 1997. Effects of the null 4A allele in near isogenic pairs of lines of wheat as expressed by starch composition and paste viscosity. In Cereals '97, Proc. 47th Australian Cereal Chem. Conf., eds A.W. Tarr, A.S. Ross and C.W. Wrigley, pp. 457-460. Royal Aust. Chem. Inst., Melbourne.
- Wootton, M., Panozzo, J.F., and Hong, S.-H. 1998. Differences in gelatinisation behaviour between starches from Australian wheat cultivars. *Starch* 50(4):147-153.
- Yun, S.H., and Quail, K. 1999. RVA pasting properties of Australian wheat starches. *Starch* 51(8-9):274-280.
- Zhu, J., O'Brien, L.O., Mares, D.J., and Shah, S.H. 1997. Quality characteristics of Chinese wheats. In: Cereals '97, Proc. 47th Australian Cereal Chem. Conf., eds A.W. Tarr, A.S. Ross and C.W. Wrigley, pp. 276-279. Royal Aust. Chem. Inst., Melbourne.