

RVA world

The technical journal of Newport Scientific

OCTOBER 2001

ISSUE 19

APPLICATION:

Dairy Product Processing in the RVA™

L. Bennett, H. Harte and G. Smithers

Food Science Australia, Werribee, Victoria 3030 Australia



The RVA is potentially useful for simulating processing of dairy products such as soft cheeses, recombined products and chilled desserts, including ice-cream. This application

describes the use of the RVA for simulating yogurt manufacturing conditions and determining the final viscosity of the fermented yogurt. Physical properties of yogurt can be affected by variability in milk

composition and processing conditions. This application can therefore be useful for quality control of ingredients including cultures, quality assurance of products, troubleshooting ingredients, or processing problems associated with yogurt making.

Stirred yoghurts were prepared from low, medium or high heat-treated skimmed milk powders (SMP), as follows: 3.0 g of SMP, 0.75 g of anhydrous milkfat and 21 g of water were mixed and warmed to 40°C before homogenisation. The quantities of ingredients should be adjusted to yield the exact

proportions of total protein, fat and moisture, as required. Exactly 15 g of the homogenised mixture was weighed into a sterile canister and pasteurised using method (a) in Table 1. The sample was inoculated by adding an aliquot (100 µl) of culture into the canister and mixing gently using the paddle. The canister and

IN THIS ISSUE...

- Stirring number at high altitude
- Tips & tricks
- Recent references


paddle were sealed with parafilm and incubated at 43°C for approximately 6 hrs or until fermentation had progressed as required. After fermentation, the yoghurt was processed using method (b) in Table 1 and the final viscosity was calculated from the mean of the final 20 readings after equilibration to 10°C. The coefficient of repeatability of the method was 8.8%. The method may be modified to more closely simulate commercial formulation and processing practices, and additional quality indices of the yogurt can be measured, as required.

The changes in viscosity during the finishing cycle demonstrate the influence of processing histories of the SMPs on the final viscosity of the yogurt (Figure 1).

TABLE 1. RVA METHODS FOR (A) YOGURT PROCESSING AND (B) TESTING FINAL VISCOSITY.

	Time	Type	Value
(a)	00:00:00	Temp.	25°C
	00:00:00	Speed	500 rpm
	00:03:00	Temp.	70°C
	00:03:00	Speed	150 rpm
	00:05:00	Temp.	70°C
	00:07:00	Temp.	95°C
	00:11:00	Temp.	95°C
	00:14:00	Temp.	43°C
	00:17:00	Temp.	43°C
	00:17:00	Speed	0 rpm
(b)	00:00:00	Temp.	43°C
	00:00:00	Speed	500 rpm
	00:00:10	Temp.	43°C
	00:00:10	Speed	150 rpm
	00:25:00	Temp.	20°C
	00:50:00	Temp.	10°C
	00:50:00	Temp.	10°C
	00:55:00	Speed	0 rpm

Figure 1. RVA pasteurisation and finishing profiles for yogurts made from low (A), medium (B) and high (C) heat SMPs.



Stirring Number Method at High Altitude

High altitude in South Africa (1500-1650 m) necessitates that ICC Method No. 161, for the determination of Stirring Number to estimate weather damage of wheat, must be modified to cook the paste at 91°C instead of the conventional 95°C so the sample does not boil.

PROFILE		
Time	Type	Value
00:00:00	Temp.	91°C
00:00:00	Speed	960 rpm
00:00:10	Speed	160 rpm

Idle temperature: 91°C ± 1°C
End of test: 3 minutes
Time between readings: 2 seconds

Figure 1. Example of SN 95°C Method in Sydney





Figure 2. Example of SN 91°C Method in Sydney



Repeatability

The method was trialed on RVA-Mini3 at the Small Grains Institute in South Africa during the 2000 grain harvest using weather damaged wheat (Table 1). A second experiment tested repeatability using sound wheat (Table 2). Results were collected in centipoise and automatically converted by the RVA to Stirring Number and Falling Number Equivalents.

Correlation


An excellent correlation with Falling Number Method was also developed during the trial:

$$\text{Lab. FN} = 35.79 + 0.1072 \times \text{cP}, R^2 = 86.8\%$$

Figure 3. Correlation with Falling Number Method

Variable	N	Mean	Median	TrMean	StDev	SE Mean
FNE Morning	20	162.86	162.80	162.96	3.70	0.83
FNE Afternoon	20	158.22	158.09	158.22	2.44	0.55
FNE Next Morning	20	158.38	158.76	158.34	3.86	0.86

Variable	N	Mean	Median	TrMean	StDev	SE Mean
Falling Number	30	265.83	268.00	265.65	14.93	2.72
FNE (SGI)	30	287.73	290.53	287.19	14.97	2.73



References:

Elliott, B., M. L. Bason, A. Leung, and L. Odendaal. Field Experience with Stirring Number: Falling Number Correlation Using a Modified ICC and AACC Method. Poster, 2nd International Wheat Quality Conference, Manhattan, KS, USA, May 2001.

Booth, R., M. Blaauw, K. Patel, and M. Bason. Detecting Sprout Damaged Wheat Using the Rapid Visco Analyser-Mini3 — A South African Experience. Poster, Pre-Harvest Sprouting Conference, Kruger National Park, South Africa, June 2001.

TIPS & TRICKS

Solving Sticky Problems in the RVA™

Starch doesn't always behave just the way you think it will in the RVA, and sometimes a little chemical assistance can be called for. Use of chemicals will usually affect viscosity results, so chemical type, amount and sequence of addition should be standardised for comparative tests.

- To moisten the sample, 2-4 g of ethanol or propanol before adding water aids the dispersion of pre-cooked starchy samples and stops the formation of lumps. Try shaking the can with an inverted stopper before starting the test to help the dispersion. If bubbles form during the test (often seen as 'jags' in the curve), delete the fast mix from the start of the RVA cycle, and remember that alcohol reduces boiling temperature, so reduce the maximum profile temperature to 90°C.¹
- To stop overflow of foaming samples from the can during the test, add 1-2 drops of anti-foaming agent such as hex-, sept- or oct- alcohols. It also helps to remove the fast mix from the profile.
- The addition of bases such as NaOH and KOH up to pH13 will help samples to gelatinise at a lower temperature. You will need to use coated RVA canisters and gloves and a fume cabinet when handling bases.
- Calcium chloride 3 Molar concentration will help samples to gelatinise at a lower temperature. This is useful for testing samples that are difficult to cook in water, such as high amylose or oxidised samples.
- To aid the dispersion of pre-cooked starchy samples and to stop the formation of lumps, moisten the sample with 2-4 g of ethanol or propanol before adding water.
- The addition of 90% aqueous dimethyl sulphoxide (DMSO) will dissolve starch without granule swelling (useful to separate swelling and dissolution effects). Use gloves and a fume cabinet when handling DMSO.²

References

- 1 Whalen, P.J., M. L. Bason, R. I. Booth, C. E. Walker, and P. J. Williams. 1997. Measurement of extrusion effects by viscosity profile using the Rapid ViscoAnalyser. *Cereal Foods World* 42(6):469-475 (Newport original 374).
- 2 Cooreman, F. L., H. van Rensburg, and J.A. Delcour, 1995. Pasting profiles and solubility of native and cross-linked corn starch in dimethylsulfoxide-water mixtures. *J. Cereal Sci.* 22:251-257 (Newport original 212).

STOP PRESS

RVA™ — The Latest Tool for Dairy Product Research

The RVA-4 is now being used for cheese sauce and other dairy product research at the University of Minnesota Department of Food Science and Nutrition (<http://fscn.che.umn.edu/>).

Assistant Professor Lloyd Metzger is something of an expert when it comes to the manufacturing, functional properties and methods of enhancing the nutritional benefits of cheese and fermented dairy products. His paper at the American Dairy Science Association Meeting, in July 2001, showed just what is possible when an RVA-4, American cheese and the University of Minnesota get together.

FOR MORE INFORMATION CONTACT:

Lloyd Metzger
Assistant Professor
University of Minnesota
Department of Food Science and Nutrition
1334 Eckles Avenue, St Paul, MN 55108-6099
tel 612-624-7786
fax 612-625-5272
e-mail lmetzger@umn.edu

RECENT REFERENCES MENTIONING THE RVA™

- Bay-Smidt, A. M., A. Blennow, M. Bojko, and B. L. Møller. 1999. The effect of phosphate and amylopectin molecular size on peak viscosity of starch pastes. *Annual Transactions of the Nordic Rheology Society*, edited by A. Friis, Vol. 7, pp. 31-38.
- Bello, A. B. 2000. Viscometry of starch in dehydrated potatoes. *Cereal Foods World* 45(10):460-463.
- Cho, K-Y., and S-T. Lim. 1998. Preparation and properties of benzyl corn starches. *Starch* 50(6):250-257.
- Dowd, M. K., M. Radosavljevic, and J. Jane. 1999. Characterization of starch recovered from wet-milled corn fiber. *Cereal Chem.* 76(1):3-5.
- Farhat, I. A., O. Tunde, and R. J. Neale. 1999. Characterisation of starches from West African yams. *J. Sci. Food Agric.* 79:2105-2112.
- Field, J. S., M.V. Swain, and N. Phan-Thien. 1996. An experimental investigation of the use of random squeezing to determine the complex modulus of viscoelastic fluids. *J. Non-Newtonian Fluid Mech.* 65:177-194.
- Grant, L.A. 1998. Effects of starch isolation, drying, and grinding techniques on its gelatinization and retrogradation properties. *Cereal Chem.* 75(5):590-594.
- Hill, S. E., and D.A. Gray, 1999. Effect of sulphite and propyl gallate or ferulic acid on the thermal depolymerisation of food polysaccharides. *J. Sci. Food Agric.* 79:471-475.
- Htoon, A. K., and A. J. Evans. 1997. The functional properties of legume starches. In *Cereals 97, Proc. 47th Australian Cereal Chem. Conf.*, edited by A. W. Tarr, A. S. Ross and C. W. Wrigley, pp. 148-150. Royal Aust. Chem. Inst., Melbourne.
- Lin, P., and Z. Czuchajowska. 1998. Role of phosphorus in viscosity, gelatinization, and retrogradation of starch. *Cereal Chem.* 75(5):705-709.
- Lizuka, K., and T. Aishima. 1999. Starch gelation process observed by FT-IR/ATR spectrometry with multivariate data analysis. *J. Food Sci.* 64(4):653-658.
- Lobeira, R., H. D. Almeida-Dominguez, and L. W. Rooney. 1998. Methods to evaluate hydration and mixing properties of nixtamalized corn flours. *Cereal Chem.* 75(4):417-420.
- McPherson, A. E., T. B. Bailey, and J. Jane. 2000. Extrusion of cross-linked hydroxypropylated corn starches. 1. Pasting properties. *Cereal Chem.* 77(3):320-325.
- Mijland, P. J. H. C., A. M. Janssen, and P. A. M. de Vres. Multivariate comparison of the Brabender Viscograph and the Rapid Visco Analyser using cross-linked starch. *Starch* 51(1):33-39.
- Moheno-Pérez, J. A., H. D. Almeida-Dominguez, and S. O. Serna-Saldivar. 1997. Effect of fiber degrading enzymes on wet milling and starch properties of different types of sorghums and maize. *Starch* 49(12):480-484.
- Nelles, E. M., J. Dewar, M. L. Bason, and J. R. N Taylor. 2000. Maize starch biphasic pasting curves. *J. Cereal Sci.* 31:287-294.
- Noda, T., Y. Takahata, T. Sato, H. Ikoma, and H. Mochida. 1996. Physicochemical properties of starches from purple and orange fleshed sweet potato roots at two levels of fertilizer. *Starch* 48(11/12): 395-399.
- Noda, T., Y. Takahata, T. Sato, H. Ikoma, and H. Mochida. 1997. Combined effects of planting and harvesting dates on starch properties of sweet potato roots. *Carbohydrate Polymers* 33:169-176.
- Perez, E. E., W. M. Breene, and Y. A. Bahnssey. 1998. Gelatinization profiles of Peruvian carrot, cocoyam and potato starches as measured with the Brabender Viscoamylograph, Rapid Visco-Analyser, and differential scanning calorimeter. *Starch* 50(1):14-16.
- Qian, J., P. Rayas-Duarte, and L. Grant. 1998. Partial characterization of Buckwheat (*Fagopyrum esulentum*) starch. *Cereal Chem.* 75(3):365-373.
- Rogers, R. and C. Sparks. 1999. The correlation of classical starch fluidity analyses with pasting curve analyses by Rapid Visco Analyser. *RVA World* 14:2-4.
- Sriburi, P., S. E. Hill, and J. R. Mitchell. 1999. Effects of L-ascorbic acid on the conversion of cassava starch. *Food Hydrocolloids* 13:177-183.
- Sriroth, K., S. Wanlapatit, K. Piyachomkwan, and C. G. Oates. 1998. Improved Cassava starch granule stability in the presence of sulphur dioxide. *Starch* 50(11-12): 466-473.
- Suhendro, E. L., C. M. McDonough, L. W. Rooney, R. D. Waniska, and S. Yetneberk. 1998. Effects of processing conditions and sorghum cultivar on alkaline-processed snacks. *Cereal Chem.* 75(2):187-193.
- Wrigley, C. W. 1999. Australian experiences in carrying research innovations to world markets. Presented at Innovations from Research in the Food Sector, European Union Conference, Trolleholm, Sweden, 11 March 1999.